



The Sonesta 6210: Engineered for optimal patient care

Positioned to be the best

- Three motors offer unparalleled positioning/access options
- Radiolucent carbon fiber tabletop for head-to-toe imaging
- 4 Double wheel rolling casters offer unencumbered manoeuvrability

The 6210 fluoroscopy procedure table offers the physician standing, seated, and supine positioning capabilities. This versatile table easily adjusts from/to standing, seated, and supine positions via a hand held control... thus creating a table that is easily accessible to patients of any age. The side-mounted cantilevered design offers the physician anterior/posterior, lateral and oblique views without obstruction. The superb c-arm access and a wide range of positioning possibilities, make the 6210 the ideal table for any fluoroscopy procedure.



SPECIFICATIONS

Length without Foot rest/leg extension:	156cm	63"
Length with Foot rest/leg extension:	196cm	78"
Width (total):	99cm	39"
Height (Min):	50cm	20"
Height (Max):	130cm	51"
Imaging area without Foot rest/leg extension (length):	156cm	63"
Imaging area with Foot rest/leg extension (length):	196cm	78"
Weight capacity:	250kg	551lbs

Standard Accessories:

- Hand control 524-6356
- 1 plastic funnel w/drainage 6325-9
- 1 Fixed arm rest 524-6021-6

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601

All Products carry the CE mark.

Images may include some optional items. Contact your representative for available configurations. Specifications subject to change without notice.