



Consumables **catalogue**

www.mmtsystems.com

FEEL THE FLOW

Urodynamics made easy with MMT






Dear Valued Customers,


MEDKONSULT medical technology s.r.o. is pleased to introduce you our disposables and accessories catalogue.

Our product list includes a complete range of disposables to provide complete urodynamic examinations in your daily practice.



Fluid filled Catheters

Abdominal Pressure Catheters



Single lumen rectal catheter



Product code	MPB4221
Tube size	12Fr (4 mm)
Tube length	300 mm
Balloon diameter	11 mm
Balloon length	55 mm
Balloon material	Latex
Qty/box	10



Double lumen rectal catheter



Product code	MBP4591
Tube size	9Fr (3 mm)
Tube length	200 mm
Balloon diameter	7 mm
Balloon length	20 mm
Balloon material	Silicone
Qty/box	10

Double lumen catheters for pressure/flow study (PFS) and cystometry




Product code	MPB5261
Tube size	6Fr (2 mm)
Tube length	400 mm
Lumens no.	2
Distance between holes	10 mm
Tip	Straight
Qty/box	10



Product code	MPB5281
Tube size	8Fr (2,67 mm)
Tube length	400 mm
Lumens no	2
Distance between holes	10 mm
Tip	Straight
Qty/box	10




Product code	MPB5481
Tube size	8Fr (2,67 mm)
Tube length	400 mm
Lumens no	2
Distance between holes	10 mm
Tip	Tiemann
Qty/box	10



Product code	MPB5661
Tube size	6Fr (2 mm)
Tube length	400 mm
Lumens no.	2
Distance between holes	10 mm
Tip	Pigtail
Qty/box	10

Triple lumen catheter for cystometry and urethral pressure profile (UPP)




Product code	MPB5804
Tube size	9Fr (3 mm)
Tube length	400 mm
Lumens no.	3
Distance between holes	50 mm
Tip	Straight
Qty/box	10

Double lumen catheter for cystometry and urethral pressure profile (UPP)



Product code	MPB5081
Tube size	8Fr (2,67 mm)
Tube length	400 mm
Lumens no	2
Distance between holes	60 mm
Tip	Straigth
Qty/box	10




Product code	MPB5891
Tube size	9Fr (3 mm)
Tube length	400 mm
Lumens no.	3
Distance between holes	10 mm, 60 mm
Tip	Straigth
Qty/box	10

Code	Name	Size (Fr)	Length (mm)
MBP4052	1L rectal catheter, PVC balloon 11x25, slit	4	1800
MPB4122	1L rectal catheter, Latex balloon, 15x75	12	300
MBP4221	1L rectal catheter, Latex balloon, 11x55	12	400
MBP4121	2L rectal catheter, Latex balloon 11x20	12	300
MBP4061	2L rectal catheter, Latex balloon, 11x25	6	400
MBP4391	2L rectal catheter, Latex balloon, 7x20	9	150
MBP4421	2L rectal catheter, Latex balloon, 11x55	12	300
MBP4441	2L rectal catheter, Latex balloon, 15x75	12	300
MBP4591	2L rectal catheter, Silicone balloon, 7x20	9	200
MBP4521	2L rectal catheter, Silicone balloon, 11x55	12	420
MBP4621	2L rectal catheter, Silicone balloon, 11x55	12	300
MBP4641	2L rectal catheter, Silicone balloon, 15x75	12	300
MBP5061	2 lumen Pves + Pura catheter; 50 mm	6	400
MBP5081	2 lumen Pves + Pura catheter; 50 mm	8	400
MBP5101	2 lumen Pves + Pura catheter; 50 mm	10	400
MPB5261	2 lumen Pves + Inf (cysto + PQ) catheter; 10 mm	6	400
MBP5281	2 lumen Pves + Inf (cysto + PQ) catheter; 10 mm	8	400
MBP5301	2 lumen Pves + Inf (cysto + PQ) catheter; 10 mm	10	400
MBP5481	2 lumen Pves + Inf (cysto + PQ) catheter; Tiemann	8	400
MBP5501	2 lumen Pves + Inf (cysto + PQ) catheter, Tiemann	10	400
MBP5661	2 lumen Pves + Inf (cysto + PQ) catheter, Pigtail	8	400
MBP5641	1 lumen Pves catheter, Pigtail	4	300
MBP5381	2 lumen Pves + UPP Radiopaque catheter; 50 mm	12	400
MBP5804	3 lumen Pves + Pura (cysto + UPP) catheter, 50 mm	9	400
MBP5802	3 lumen Pves + Pura (cysto + UPP) catheter, Tiemann	9	400
MBP5891	3 lumen Pves + Pura (cysto + UPP) catheter; 10/60 mm	9	400

- Manufacturer: Bioengineering Laboratories s.r.l (BEL), Italy
- • Highlighted items are in stock, other items with 4-6 weeks delivery

Abdominal Pressure Catheters



Single lumen rectal catheter


←

Product code	MA001145
Tube size	4,5Fr
Tube length	2000 mm
Balloon type	Slit
Balloon length	20 mm
Balloon material	PVC
Qty/box	25

Double lumen rectal catheter


→

Product code	MA002110
Tube size	10Fr
Tube length	385 mm
Extension lines	4,5Fr x 1300 mm; 6Fr x 150 mm
Balloon material	Poly-isoprene
Qty/box	10



Code	Name	Size (Fr)	Length (mm)
MA001145	1L rectal catheter, PVC balloon, slit	4,5	2000
MA002110	2L rectal catheter, Poly-isoprene balloon	10	385
MA011145	1L rectal catheter, Poly-isoprene balloon	4,5	2000
MA051145	1L rectal catheter with introducer, PVC balloon	4,5	2000
MA061145	1L rectal catheter with introducer, Poly-isoprene balloon	4,5	2000
MA052110	2L rectal catheter, PVC balloon	10	2000
MA101008	Bladder filling catheter	8	385
MA101145	Bladder pressure catheter	4,5	2000

- Manufacturer: PB Medical Ltd, UK
- ● Highlighted items are in stock, other items with 6 - 8 weeks delivery



Extension set 2 Extension set 3



Product code	MBP1282 / MBP1292
Internal diameter	1,5 mm
Length	1,5 m
Material	PVC
Type	Disposable
Connector color	Blue and Red / Blue, red and Clear

Pressure transducer



Product code MA9201-1


Disposable pressure transducer with 3-way stopcock and 1-way stopcock.



Urodynamic set for pump



Product code	MU9109
Length	1845 mm
Middle part material	Silicone
Type	Disposable
Qty/box	25



EMG electrodes



Product code	9106
Qty/package	12
Type	Disposable



Y adapter for profilometry



Product code	9109
Material	Transparent PVC tube
Length	Silicone
Type	Disposable
Qty/box	30



MEDKONSULT medical technology s.r.o.
Pasteurova 15, 779 00 Olomouc, Czech Republic

T +420 581 113 030
E info@mmtsystems.com

www.mmtsystems.com