

Product catalogue 2020



© 2020 Sooma Oy, Helsinki, Finland. All rights reserved. Reproduction in any part without the written permission of Sooma Oy is prohibited.

Sooma reserves the right to make changes in any specifications shown herein at any time without notice or obligation. All specifications are for informational purposes only.

Sooma and Sooma tDCS are registered trademarks of Sooma Oy. Pending international patents cover the equipment and software that comprise the Sooma tDCS™ system.

Table of contents	3
Sooma Depression Therapy	4
Sooma Pain Therapy	5
ComfoCaps	6
ComfoTodes and ComfoLeads	7
ComfoPads	8
Support products and service	9

CE
0598



Sooma Depression Therapy is a solution for bifrontal transcranial direct current stimulation. It includes a stimulation system and a remote patient monitoring software.

The Home Kit includes all you need for administering Sooma Depression Therapy efficiently both in the clinic or at the patients' home.

The stimulator unit comes with 2 mA current and 30-minute stimulation duration as default parameters. This current is applied through the DLPFC using proprietary Sooma Comfort accessories.

The solution is Class IIa medical device indicated for treatment of major depressive disorder. It is built to the highest industry standards in Finland.

Order code	Item	Price
DCS2H	<p>Sooma Depression Therapy Home Kit</p> <ul style="list-style-type: none"> o Sooma tDCS™ stimulator configured for Depression Therapy o Access to Sooma Software Suite o ComfoTrode cups (CE25) o ComfoLeads cables (CSNP) o ComfoCap, Medium. Montage 2: Bifrontal stimulation (CCM2) o 20 pairs of ComfoPads (CP25) o Neck Pouch for Sooma tDCS™ (NECK) o Sooma tDCS™ Mobile Case (CAMO) o 12-month warranty 	

CE
0598



Sooma Pain Therapy is a solution for transcranial direct current stimulation of the motor cortex. It includes a stimulation system and a remote patient monitoring software.

The Home Kit includes all you need for administering Sooma Pain Therapy efficiently both in the clinic or at the patients' home.

The stimulator unit comes with 2 mA current and 20-minute stimulation duration as default parameters. This current is applied to the motor cortex using proprietary Sooma Comfort accessories.

The solution is Class IIa medical device indicated for treatment of Fibromyalgia and Chronic Neuropathic Pain. It is built to the highest industry standards in Finland.

Order code	Item	Price
DCS3H	Sooma Pain Therapy Home Kit <ul style="list-style-type: none"> o Sooma tDCS™ stimulator configured for Pain Therapy o Access to Sooma Software Suite o ComfoTrode cups (CE25) o ComfoLeads cables (CSNP) o ComfoCap, Medium. Montage 3: Motor cortex stimulation (CCM3) o 20 pairs of ComfoPads (CP25) o Neck Pouch for Sooma tDCS™ (NECK) o Sooma tDCS™ Mobile Case (CAMO) o 12-month warranty 	



Sooma ComfoCap series provides an intuitive and accurate method for electrode placement. It offers an industry leading reproducibility in positioning allowing every stimulation session to maximize the therapeutic efficacy.

ComfoCaps incorporate a hydrophobic fabric mixture that protects from electrode bridging. This allows you to know that you are always reaching the stimulation target.

There are standard variants of Sooma ComfoCap for bifrontal stimulation (Montage 2) and motor cortex stimulation (Montage 3). The Montage 3 ComfoCap includes a capability to change the stimulation polarity. Additionally, ComfoCaps can be ordered without default placement allowing you to create a patient-specific montage.

The lifetime for ComfoCaps is 3 months or 100 sessions.

Order code	Item	Price
CCS2 CCM2 CCL2	ComfoCap, Montage 2 for bifrontal (F3 + F4) stimulation <ul style="list-style-type: none"> ○ Small, Head circumference: 50-54 cm ○ Medium, 55-58 cm ○ Large, 59- cm 	
CCS3 CCM3 CCL3	ComfoCap, Montage 3 for motor cortex (C3 + Fp2 / C4 + Fp1) stimulation <ul style="list-style-type: none"> ○ Small, Head circumference: 50-54 cm ○ Medium, 55-58 cm ○ Large, 59- cm 	
CCS0 CCM0 CCL0	ComfoCap, Montage 0 for custom stimulation placement <ul style="list-style-type: none"> ○ Small, Head circumference: 50-54 cm ○ Medium, 55-58 cm ○ Large, 59- cm 	

Sooma ComfoTrodes are high-quality stimulation electrodes for transcranial direct current stimulation.



Together with ComfoPads, these electrodes offer unrivalled dose precision: The ComfoTrode design ensures that the stimulation is not directed elsewhere by excess electrolyte, but goes only to the intended stimulation area.

The lifetime for ComfoTrodes is 3 months or 100 sessions.

Order code	Item	Price
CE25	ComfoTrode stimulation electrode cup, round 25cm2	



ComfoLeads are touch protected cables for connecting ComfoTrodes to your stimulation unit. The 1.5 mm cables and Snap-On connectors are color-coded.

The lifetime for ComfoLeads is 12 months.

Order code	Item	Price
CSNP	ComfoLeads cables, 150mm color-coded red - black	

CE



Meet ComfoPad – The conductive surface with industry leading comfort.

ComfoPads expand rapidly from a dry state when soaked in electrolyte to create a soft and uniform stimulation surface. One pad can absorb up to 15ml of electrolyte and release it only when pressed against the skin surface. This creates up to 40% better contact quality when compared to standard stimulation sponges.

The enhanced contact quality offered by ComfoPad translates directly into weaker stimulation sensation. Consequently, stimulation with higher intensities becomes possible and reliability of placebo-control in crossover setting greatly improves.

Finally, use of ComfoPad saves several minutes of preparation time each session.

ComfoPads are single-use only.

Order code	Item	Price
CP25	ComfoPad, swelling stimulation pad - For 25cm ² electrode, 20 pairs	

Electrolyte



Optimal electrolyte solution for tDCS: Sodium chloride in 30ml bottles suitable for single session.

Features

- o Five 30ml bottles of sterile medical grade NaCl
- o Only Na⁺ and Cl⁻ ions
- o NaCl content 0.9%

Mobile Case



Mobile case for Sooma tDCS™ device and accessories. Fits everything needed for one homebased therapy program.

Features

- o Enables safe and easy transportation of the device and accessories
- o Keeps the device and accessories safe in-between sessions

Order code	Item	Price
NACL	Stimulation electrolyte, NaCl 0.9%. For 5 stimulation sessions (5x 30 ml)	
CAMO	Sooma tDCS mobile case	
SERV	Service or configuration work. Sooma tDCS system offers the following configuration possibilities: <ul style="list-style-type: none"> o Output current: Default 2.0 mA. Range 0-3.0 mA (0.1mA steps) o Stimulation duration: Default 30 mins. Range 0-60 (30 second steps) o Mode: Active, Sham o Treatment lock (15 sessions): OFF, ON o Language: EN, DE, FR, IT, SP, PT, PL, SE or FI 	
WA36	36-month warranty extension	
WA60	60-month warranty extension	



Sooma Oy
Kuortaneenkatu 2
00510 Helsinki
FINLAND

www.soomamedical.com



+358 10 328 9811



sales@soomamedical.com